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### \* *1M Tris-1 Liter*

Dissolve 121.14 g Tris (hydroxymethyl) aminomethane (tris, MW= 121.14 in 800 ml H<sub>2</sub>O). Adjust PH

to the desired value by adding concentrated HCL as following table;

1M Tris with PH 7.4 need approximately 70 ml of concentrated HCL

1M Tris with PH 7.6 need approximately 60 ml of concentrated HCL

1M Tris with PH 8.0 need approximately 42 ml of concentrated HCL

Adjust volume to 1 liter with H<sub>2</sub>O. Sterilize by autoclaving and store at room temperature

A series of horizontal dotted lines spanning the width of the page, providing a template for handwritten notes.

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Adjust to PH 8.0 with NaOH (20 g NaOH pellets) and adjust volume to 1 liter with  $\text{H}_2\text{O}$ .

Divide into aliquots and sterilize by autoclaving.

The disodium salt of EDTA will only solute when the PH of the solution is adjusted to 8.0 by the addition of NaOH.

A series of horizontal dotted lines spanning the width of the page, providing a template for handwritten notes.

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**\* 1 X TE -1 litter**

Add 10 ml 1 M Tris (PH 8.0, 7.6 or 7.4) and 2 ml 0.5 M Na<sub>2</sub> EDTA (PH 8.0) to 800 ml H<sub>2</sub>O. Mix and adjust volume to 1 liter with H<sub>2</sub>O. Sterilize by autoclaving .Store at room temperature.

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**\* 3M NaOAc**

Dissolve 408.24 g sodium acetate-3H<sub>2</sub>O (NaOAc-3H<sub>2</sub>O, MW-136.08) in 800 ml H<sub>2</sub>O. Adjust to PH 5.2 with glacial acetic acid or to PH 7.0 with diluted acetic acid.  
Adjust volume to 1 liter with H<sub>2</sub>O. Sterilize by autoclaving and store at room temperature.

A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.

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**\* 10 X TBE buffer**

Dissolve 108 g Tris and 55 g Boric acid in 900 H<sub>2</sub>O. Add 40 ml 0.5M Na<sub>2</sub>EDTA (PH 8.0) and adjust volume to 1 liter with H<sub>2</sub>O. Store at room Temperature.

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**\* 50 X TAE**

Dissolve 242 g tris in 500 ml H<sub>2</sub>O. Add 100 ml 0.5 M Na<sub>2</sub>EDTA (PH 8.0) and 57.1 ml glacial acetic acid. Adjust volume to 1 liter with H<sub>2</sub>O. Store at room temperature.

A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.

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### \* *Ethidium bromide (500µg/ml)*

Dissolve 50 mg ethidium bromide in 100 ml distilled water. Store at 4°C in a dark bottle for no longer than one year. Dilute stock 1:1000 to achieve working concentration of 0.5 µg / ml.

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### \* *Loading Buffer*

Dissolve 0.05 gram of bromophenol blue , 0.05 gram xylene cyanol , 0.186 gram Na<sub>2</sub> EDTA in 5 ml glycerol and bring the volume to 10 ml with TBE buffer.

A series of horizontal dotted lines spanning the width of the page, providing a template for handwritten notes.

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### \* *Agarose Gel Preparation (0.8%)*

Weigh out 0.8 gram agarose powder and transfer to 200 ML bottle. Then add 100 ml of TE or TBE (1X) buffer. Dissolve the agarose by heating in a microwave. Swirl to ensure that all agarose is completely dissolved. Once the temperature of agarose reached 65<sup>0</sup>C, add 100 µl of ethidium bromide. Swirl and pour the gel immediately in the tray. Let it solidify.

A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.

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**\* 20 X SSC solution**

Dissolve 175.3 g of NaCl ( NaCl M.W= 58.44) and 88.2 gram of Na-citrate (Na-citrate M.W=294.1) in 1000 ml of D.W. Bring the PH into 7.0 using 5M NaCl. For 2 X SSC, add 1 volume of 20 X SSC to 9 volumes of D.W. For 5X SSC, add 1 volume of 20 X SSC to 3 volume of D.W.

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**\* 1 M HCL**

Prepare 1 M HCL by adding 85.5 ml concentrated HCL (37%, Merck .P.a) to 914.5 ml water. Store at room temp for no longer than 6 months. Prepare 0.25 M HCL directly for use by adding 125 ml 1 M HCL in 375 ml D.W

A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.

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**\* 0.4 M NaOH**

Prepare 4 M NaOH by dissolving 160 g NaOH in 1000 ml of D.W. Store at room temperature for no longer than 3 months. Prepare 0.4 M NaOH directly for use by adding 100 ml 4M NaOH in 900 ml D.W.

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#### \* *Primary wash buffer*

Dissolve 360 gram urea, 4 gram SDS, and 25 ml of 20 X SSC in to 1000 ml D.W. Do not autoclave. Either make just before use or sterilize by filtrating through 0.2 mm filter and store at room temp for no longer than one year.

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
### \* *Estimation of DNA Concentration*

DNA concentration can be performed by staining DNA with ethidium bromide, after short electrophoresis or by using spectrophotometer.

The quality and the approximate concentration of extracted DNA were estimated before digestion by staining of DNA with ethidium bromide, after short electrophoresis (15 minutes at 85 V) on 0.8% agarose gel. (Fig 1)

A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.


Fig 1: DNA of Mycobacterium

A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.

### **\* Digestion of DNA**

DNA of mycobacteria digested using PVUII restriction endonuclease that cleaves IS6110 once. Always prepare the mix in the following order: water, restriction buffer, DNA and enzyme. (Fig 2)

Add 1  $\mu$ l of enzyme, 3  $\mu$ l of enzyme buffer, 5  $\mu$ l of DNA, and bring the volume to 20-25  $\mu$ l with D.W (depending on DNA concentration). The Mix incubated in 37<sup>0</sup>C for 18-24 hours.

Load samples (about 2  $\mu$ l) on 1% agarose in 1x TBE buffer, in order to give an idea about digestion of the samples. (Electrophoresis at 100 voltages for 20 minutes). (Fig 3)

If the DNA were digested completely, then prepare a big gel and load the whole digested samples in to each wells and run the gel overnight with voltage 25-30. (Fig 4)

Always include H<sub>37</sub>R<sub>V</sub> in gel, because it serves as an external marker for computer analysis. (Fig 3)

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

A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.

**Figure 2.1** DNA can be cleaved by restriction enzymes into fragments that can be separated by gel electrophoresis.

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
*DNA digested after short electrophoresis*

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
*DNA digested after long electrophoresis*

A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.

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### **\* *Depurination step* \***

Place the gel with its tray for 10 min in 500 ml 0.25 M HCL. Rinse the gel with distilled water. The whole step repeats once more.

A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.

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**\* Denaturation step**

Place for 15 min in 500 ml (0.4 M NaOH) in shaker. Rinse the gel with distilled water.  
Fill the tray with 0.4 M NaOH and put in shaker for 20 minutes. Rinse the gel with distilled water. At this stage the gel is ready for transfer.

A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.

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### \* *Blotting by capillary method*

Place a gel-holder upside down in a plastic tray. Fill the dish with Alkaline buffer (0.4N NaOH with 1M NaCL). Cut the Whatman filter paper exactly the same size of the gel. Place it on gel-holder in a way that the side of the paper floated into plastic tray that has been filed with alkaline buffer. Place the gel on gel-holder and cover it with wet nylon membrane (Wet the paper with D.W for 5 minutes and then soak in 0.4 N NaOH). Cover all 4 sides of the filter with parafilm. Place 6 thin Whatman papers followed by 6 thick Whatman papers on the filter. To improve the capillary transfer, the gel was compressed by placing 4 to 5 kg weight on the top of the blotting tower. The complete transfer was occurred in 18-24 hours. (Fig 5)

A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.


Fig 5: Blotting by Capillary Method

A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.

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***\* Treatment of the membrane after blotting by capillary***

After the completion of transfer, the membrane washed by Tris- solution for 15 min. Then, it wrapped up by cellophane plastic sheet and fixed using UV trans- illuminator for 5 min.

A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.

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A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.

**تھیہ DNA پروب با PCR**  
**Preparation of DNA probe by PCR**

dNTP: 1  $\mu$ l (0.2 mM)  
 INS 1: 1  $\mu$ l (20 pmol)  
 INS 2: 1  $\mu$ l (20 pmol)  
 DNA Template: 2.5  $\mu$ l (50-100)

Mgcl2: 1.5  $\mu$ l (1.5 mM)  
 PCR Buffer: 5  $\mu$ l (1 $\times$ )  
 Enzyme (taq): 0.25  $\mu$ l (1.25 U)

**\* PCR Cycle:**

<b>Initial denaturation 6 min 94<math>^{\circ}</math>c</b>			
30 Cycles	Step1	1 min	94 $^{\circ}$ C
	Step 2	1 min	61.5 $^{\circ}$ C
	Step 3	1 min	72 $^{\circ}$ C
<b>Final extension 10 min 72 <math>^{\circ}</math>C</b>			


A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.

The PCR products can be examined by electrophoresis and staining with ethidium bromide. The PCR fragments should be 245 bp in length using primer INS-1 and INS-2.

***\* Purification of the PCR product:***

For this step use the Quick Spin Sephadex G50 column (as instructed by manufacturer) or any other PCR-product purification kit.

A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.

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### \* *Detection & Labeling*

The system involves direct labeling of DNA probe with the enzyme horseradish peroxidase (HRP). This is achieved by completely denaturing the probe so that it is single-stranded and therefore negatively charged. HRP, which has been complexes with a positively charged polymer, is added and it forms a loose attachment to the DNA by charge attraction.

Addition of glutaraldehyde causes the formation of chemical cross-links so that probe is covalently bound to the enzyme. Once labeled, the probe is used in hybridization with target DNA immobilized on a membrane. HRP catalyses the oxidation of the substrate luminal which, in the presence of a chemical enhancer, results in a large and sustained light emission, that is readily detected on film.

A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.

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\* *Hybridization and Detection include four stages (using ECL detection system):*

- Labeling of the probe
- Hybridization
- Washing the membrane N+ after hybridization
- Detection

A series of horizontal dotted lines spanning the width of the page, providing a template for handwritten notes.

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Denature the probe for 5 min by placing it in a boiling water bath. Chill on ice immediately for one min.

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Spin the probe-mixture a few seconds in a micro-centrifuge and add an equal volume DNA labeling reagent. Mix well by carefully pipetting.

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Add the same volume glutaraldehyde. Mix well by carefully pipetting.

A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.

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Prewarm the primary wash buffer to 42°C.

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Place the filter in a clean plastic box.

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Add 500 ml secondary wash buffer to the filter and wash the filter for 5 min at room temp on a shaking platform.  
Repeat this stage twice.

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Make sure that all equipment is available because the work in the darkroom should not be interrupted.

A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.

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Close the cassette for exactly one min.

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
Place the film in the radiology machine.

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A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.

A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.

## References:

- \* Van Soolingen D, De Haas PE, Hermans PW, Van Embden JD. DNA fingerprinting of *Mycobacterium tuberculosis*. *Methods Enzymol.* 1994; 235:196-205. Review.
- \* Van Soolingen D, Hermans PW, de Haas PE, van Embden J. Insertion element IS1081-associated restriction fragment length polymorphisms in *Mycobacterium tuberculosis* complex species: a reliable tool for recognizing *Mycobacterium bovis* BCG. *J Clin Microbiol.* 1992 Jul; 30(7):1772-7.
- \* Van Embden JD, Van Soolingen D, Small PM, Hermans PW. Genetic markers for the epidemiology of tuberculosis. *Res Microbiol.* 1992 May; 143(4):385-91. Review. No abstract available.
- \* Farnia P, Masjedi MR, Varahram M, et al. The recent-transmission of *Mycobacterium tuberculosis* strains among Iranian and Afghan relapse cases: a DNA-fingerprinting using RFLP and spoligotyping. *BMC Infect Dis.* 2008 Aug 6; 8:109
- \* Ramazanzadeh R, Farnia P, Amirmozafari N, et al. Comparison between molecular epidemiology, geographical regions and drug resistance in *Mycobacterium tuberculosis* strains isolated from Iranian and Afghan patients. *Chemotherapy.* 2006; 52(6):316-20. Epub 2006 Sep 27
- \* Farnia P, Masjedi MR, Mirsaeidi M, et al. Prevalence of Haarlem I and Beijing types of *Mycobacterium tuberculosis* strains in Iranian and Afghan MDR-TB patients. *J Infect.* 2006 Nov; 53(5):331-6. Epub 2006 Feb 1
- \* Farnia P, Mohammadi F, Masjedi MR, et al. Evaluation of tuberculosis transmission in Tehran: using RFLP and spoligotyping methods. *J Infect.* 2004 Aug; 49(2):94-101
- \* Doustdar F, Khosravi AD, Farnia P, et al. Molecular analysis of isoniazid resistance in different genotypes of *Mycobacterium tuberculosis* isolates from Iran.
- \* Van Soolingen D, Hermans PW, De Haas PE, Soll DR, Van Embden JD. Occurrence and stability of insertion sequences in *Mycobacterium tuberculosis* complex strains: evaluation of an insertion sequence-dependent DNA polymorphism as a tool in the epidemiology of tuberculosis. *J Clin Microbiol.* 1991 Nov; 29(11):2578-86

A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.

Shaheed Beheshti University of  
Medical sciences and Health Services

Collaborating Center for Tuberculosis  
East Mediterranean Region

Mycobacteriology  
Research Center

National Research Institute  
Tuberculosis and Lung Disease

Collaborating Center for Tuberculosis  
Laboratory IUATLD

Address: Mycobacteriology Research Centre, Iranian National Reference TB Laboratory, National Research Institute of Tuberculosis & Lung Disease (NRITLD), Shaheed Beheshti University of Medical Sciences, Shaheed Bahonar Ave, Daarabad, P.C: 19569/44413, P.O Box: 19575/154, Tehran, Iran.

Tel: +98(21) 27122571

Fax: +98(21) 26109505

E-mail: [mrc @ nritld.ac.ir](mailto:mrc@nritld.ac.ir), [pfarnia @ nritld.ac.ir](mailto:pfarnia@nritld.ac.ir), Website: [www.mrc.ac.ir](http://www.mrc.ac.ir)

A series of horizontal dotted lines spanning the width of the page, intended for handwritten notes.